

# Experiment FB-1: Osmoregulation

## Background

Sodium is the predominant cation in the extracellular fluid of multicellular animals. The high level of sodium ions in seawater results in minimal osmotic stress in many marine organisms. In these animals, therefore, minimal ionic and osmotic regulation is required. Furthermore, the large volume of water in the ocean ensures minimal fluctuations of the osmotic environment.

Some marine organisms do not spend all of their time in the ocean. Some live in the intertidal region where they experience periods of drought, while others may live in tidal pools. In the latter case, the relatively small volume of water in the tidal pool may result in fluctuations in the osmotic environment. The level of sodium in the pool may increase when the sun evaporates water or decrease when freshwater is added through rain or a river. Thus, animals that are trapped in a tidal pool must be able to adapt to short-term osmotic stress and survive for about 12 hours until the next tide.

In this laboratory, you will use a series of dilutions of seawater (with deionized water) to measure the effects of solute concentration on the movement of water into or out of a polychaete worm. You will place a worm in each solution and then measure its weight change every 10 minutes for one hour.

## Equipment Required

- PC Computer
- IWX/214 data acquisition unit
- USB cable
- IWX/214 power supply
- FT-104 Force Transducer
- Ring stand and clamp
- Basket and pennies (each weighs about 3 g)
- 5 x 250 ml beakers
- 2 x 100 ml graduated cylinders.
- Forceps
- Artificial Seawater (see Appendix

## IWX/214 Setup

- 1 Place the IWX/214 on the bench, close to the computer.
- 2 Check Figure T-1-1 in the Tutorial chapter for the location of the USB port and the power socket on the IWX/214.
- 3 Check Figure T-1-2 in the Tutorial chapter for a picture of the IWX/214 power supply.
- 4 Use the USB cable to connect the computer to the USB port on the rear panel of the IWX/214.

- 5 Plug the power supply for the IWX/214 into the electrical outlet. Insert the plug on the end of the power supply cable into the labeled socket on the rear of the IWX/214. Use the power switch to turn on the unit. Confirm that the red power light is on.

## Start the Software

- 1 Click on the LabScribe shortcut on the computer's desktop to open the program. If a shortcut is not available, click on the Windows Start menu, move the cursor to **All Programs** and then to the listing for **iWorx**. Select **LabScribe** from the **iWorx submenu**. The LabScribe Main window will appear as the program opens.
- 2 On the **Main window**, pull down the **Settings menu** and select **Load Group**.
- 3 Locate the folder that contains the settings group, **IPLMv4.iwxgrp**. Select this group and click **Open**.
- 4 Pull down the **Settings menu** again. Select the **Osmoregulation-LS2** settings file.
- 5 After a short time, LabScribe will appear on the computer screen as configured by the **Osmoregulation-LS2** settings.
- 6 For your information, the settings used to configure the recording channels in the LabScribe software and IWX/214 for experiment are listed in Table FB-1-1 on page FB-1-1.

**Table FB-1-1: Settings on the Channel Window of the Preferences Dialog that Configure the iWorx System for Experiment FB-1.**

Parameter	Units/Title	Setting	Mode/Function
Acquisition Mode		Chart	
Start		User	
Stop		User	
Display Time	Sec	10	
Speed	Samples/Sec	200	
Channel A3	Weight	√	Record

## FT-104 Setup

- 1 Locate the FT-104 force transducer in the IWX/214 kit (Figure FB-1-1 on page FB-1-1).



Figure FB-1-1: The FT-104 force transducer.

- 2 Plug the DIN8 connector of the FT-104 into the Channel 3 input (Figure FB-1-2 on page FB-1-2).



Figure FB-1-2: The FT-104 force transducer connected to the IWX/214.

- 3 Attach the transducer to a ring stand using a 90° clamp, so that the transducer is horizontal (Figure FB-1-3 on page FB-1-2).
- 4 Attach a weight pan to end of the transducer arm.




Figure FB-1-3: Equipment setup to record weight placed in a pan.

## Preparation of Solutions

- 1 Place or write a label on each 250 ml beaker: 100%, 90%, 80%, 70% and 50%.

**Table FB-1-2: Composition of Solutions for Experiment**

Labeled Beaker	Artificial Seawater (ml)	Deionized Water (ml)
100%	200	0
90%	180	20
80%	160	40
70%	140	60
60%	120	80

- 2 Make up the seawater solutions (Table FB-1-2 on page FB-1-2) and place in the appropriately labeled beakers.

## Calibration of FT-104 Force Transducer

**Aim:** To calibrate the force transducer used to measure weight gain or loss of the animal.

### Procedure

- 1 Make sure that the IWX/214 is turned on and the FT-104 force transducer is connected to the DIN8 input for ten minutes before the calibration is performed.
- 2 Type **No Weight** in the **Mark** box to the right of the **Mark** button. Click the **Record** button, and press the **Enter** key to attach a comment to the recording. Record for ten seconds with no weight hanging from the arm or hook of the transducer.
- 3 Count the number of pennies you have and multiply their number by three (the weight of a penny in grams).
- 4 Type the weight of the pennies in the **Mark** box. Place the pennies on the weight pan and press the **Enter** key on the key board, simultaneously. Click the **AutoScale** button next to the channel title area. Record for ten more seconds (Figure FB-1-4 on page FB-1-2).




Figure FB-1-4: Recording of different known weights used to calibrate FT-104 force transducer.

- 5 Click **Stop** to halt the recording.
- 6 Select **Save As** in the **File** menu, type a name for the file. Choose a destination on the computer in which to save the file, like your lab group folder). Designate the file type as **\*.iwxdata**. Click on the **Save** button to save the data file.
- 7 Remove the pennies from the weight pan.

### Units Conversion

- 1 Scroll to the beginning of data when no weight was attached to the FT-104 force transducer.
- 2 Use the **Display Time** icons on the **LabScribe** toolbar (Figure FB-1-5 on page FB-1-3) to adjust the **Display Time** of the **Main window** to show the complete calibration data on the **Main window**. The required data can also be selected by:

- Placing the cursors on either side of data required.
  - Clicking the **Zoom between Cursors** button on the **LabScribe** toolbar to expand the complete calibration data to the width of the **Main window**.
- 3 Click the **Double Cursor icon** (Figure FB-1-5 on page FB-1-3) so that two blue cursors appear on the **Main window**. Place one cursor on the flat section of data collected when no weight was attached to the FT-104, and the second cursor on the flat section of data collected when the pennies were attached to the transducer.




Figure FB-1-5: The **LabScribe** toolbar.

- 4 To convert the output of the FT-104 force transducer from a voltage to the grams of force:
- Click on the arrow next to the title of the **Weight channel** to open the **channel menu**.
  - Select **Units** from the **channel menu** and **Simple** from the **Units submenu**.




Figure FB-1-6: The **Simple Units Conversion** dialog box window with the voltages at the cursors set to equal the weight used in calibration.

- 5 The **Simple Units Calibration** window will appear (Figure FB-1-6 on page FB-1-3). On this window:
- Select **2 point calibration** from the pull-down menu in the upper-left corner of the window.
  - Put a check mark in the box next to **Apply units to all blocks**.
  - Notice that the voltages from the positions of the cursors are automatically entered into the value equations.
  - Enter **zero** in the corresponding box to the right of the voltage recorded when no weight was attached to the transducer. Enter **the weight of the pennies** in the box to the right of the corresponding voltage recorded when the weight of the pennies was hung on the hook of the transducer.

- Enter the name of the units, **grams**, in box below the weights. Click on the **OK button** in the lower right corner of the window to activate the units conversion.

## Practice Weighing Technique

Aim: To develop a consistent technique of weighing a worm.

### Procedure

- 1 Use forceps to remove a worm from seawater. Blot the worm with paper towels to remove excess water.
- 2 Type **No Worm** in the **Mark box** to the right of the **Mark button**.
- 3 Click the **Record button** to record a baseline of ten seconds while only the weight pan is attached to the transducer. Continue recording.
- 4 Type **Worm** in the **Mark box**. Place the worm on the weight pan attached to the transducer.
- 5 Continue to record for ten seconds after the worm was placed on the weight pan. Click **Stop** to halt the recording.
- 6 Replace the worm in the seawater.
- 7 Repeat Steps 1 through 6 for two other worms.

### Data Analysis

- 1 Scroll through the data file to the weighing of the first worm.
- 2 Use the **Display Time** icons to adjust the **Display Time** of the **Main window** so that the output of the transducer before and after the worm was placed on the weight pan is displayed on the **Main window**. This section of data can also be selected by:

- Placing the cursors on the two levels of the transducer output before and after the worm was placed on the weight pan, and
- Clicking the **Zoom between Cursors** button on the **LabScribe** toolbar (Figure FB-1-5 on page FB-1-3) to expand or contract this section of recording to the width of the **Main window**.

- 3 Click on the **Analysis window** icon in the toolbar (Figure FB-1-5 on page FB-1-3) or select **Analysis** from the **Windows menu** to transfer the data displayed in the **Main window** to the **Analysis window**.
- 4 Look at the **Function Table** that is above the **Weight Channel** in the **Analysis** window. The mathematical functions, **V2-V1** and **T2-T1**, should appear in this table. The values for these parameters are displayed in the table across the top margin of the **Weight channel**.
- 5 Maximize the height of the trace on the **Weight Channel** by clicking on the arrow to the left of the channel's title to open the **channel menu**. Select **Scale** from the menu and **AutoScale** from the **Scale submenu** to increase the height of the data on that channel.
- 6 Place the cursors on the recording before and after the addition of the worm to the weight pan. The value for **V2-V1** on the **Weight channel** is the weight of the worm.

- 7 The weight of the worms can be recorded in the on-line notebook of LabScribe by typing their names and values directly into the **Journal**.
- 8 The functions in the **channel pull-down menus** of the **Analysis window** can also be used to enter the name and value for **V2-V1** from the recording to the **Journal**. To use these functions:
  - Place the cursors at the locations used to measure the weight of the worm.
  - Transfer the name of the mathematical function used to determine the weight to the **Journal** using the **Add Title to Journal** function in the **Weight Channel pull-down menu**.
  - Transfer the weight to the **Journal** using the **Add Ch. Data to Journal** function in the **Weight Channel pull-down menu**.
- 9 Repeat Steps 2 through 8 to find the weights of the other two worms weighed in this exercise. Record the weights in the **Journal** and on Figure FB-1-3 on page FB-1-4.
- 10 Select **Save** in the **File menu**.

**Table FB-1-3: Weight of Worms Recorded During Weighing Practice.**

	Weights (g)		
	Worm 1	Worm 2	Worm 3
Practice			

## Osmoregulation

**Aim:** To measure changes in the weights of five (5) worms placed in different osmotic environments using a transducer.

### Procedure

- 1 Use forceps to remove a worm from seawater. Blot the worm with paper towels to remove excess water.
- 2 Type **No Worm** in the **Mark box** to the right of the **Mark button**.
- 3 Click the **Record button** to record a baseline of ten seconds while only the weight pan is attached to the transducer. Continue recording.
- 4 Type **Worm 1** in the **Mark box**. Place the worm on the weight pan attached to the transducer. Press the **Enter key** on the keyboard to mark the recording.
- 5 Continue to record for ten seconds after the worm was placed on the weight pan. Click **Stop** to halt the recording. Note the time when the worm was weighed.

- 6 Record the weight of the worm in the **Journal** and on Table FB-1-4 on page FB-1-4.
- 7 Place the worm in the 100% seawater solution.
- 8 Repeat Steps 1 through 6 for each of the four other worms used in this exercise. Each worm goes into a different solution: 90, 80, 70, or 60% seawater.
- 9 Every ten minutes, remove each worm from its solution. Blot the worm, weigh it, and return it to the same solution.
- 10 Record the weight of the worm in the **Journal** and on Table FB-1-4 on page FB-1-4.
- 11 Weigh all the worms until you have seven weight values for each worm.

**Table FB-1-4: Weights of Worms Placed in Solutions of Different Osmotic Strengths Over Time.**

Solutions	Weights (g)				
	100%	90%	80%	70%	60%
Time (min)	Worm 1	Worm 2	Worm 3	Worm 4	Worm 5
0					
10					
20					
30					
40					
50					
60					

### Data Analysis

Graph the weight of each of the five worms in a different solution as a function of time.

### Questions

- 1 Does the weight of the worm in 100% seawater change? Is the weighing of the worms accurate?
- 2 Which worm gained weight at the fastest rate? If weight gain indicates water intake, explain the results in terms of concentration gradients?
- 3 Do any of the worms stop gaining weight towards the end of the experiment? How do you explain this observation?

## Appendix

**Table FB-1-5: Recipe for Artificial Seawater**

Concentration (mMolar)	Salt	Grams/Liter DI H <sub>2</sub> O
425.0	Sodium Chloride	24.6
9.0	Potassium Chloride	0.67
25.5	Magnesium Sulfate*7H <sub>2</sub> O	6.29
23.0	Magnesium Chloride*6H <sub>2</sub> O	4.66
2.0	Sodium Bicarbonate	0.18
9.3	Calcium Chloride*2H <sub>2</sub> O	1.36

Dissolve all the salts, except CaCl<sub>2</sub>\*2H<sub>2</sub>O, in 750 ml of deionized (DI) H<sub>2</sub>O. Dissolve CaCl<sub>2</sub>\*2H<sub>2</sub>O in 200 ml DI H<sub>2</sub>O before adding it to the solution containing the other salts to prevent precipitation of calcium carbonate. Adjust the pH of the artificial seawater to 8.0 and add DI H<sub>2</sub>O to bring the final volume to 1 Liter.