

# iWorx Physiology Lab Experiment

Experiment HE-9

## Resting, Active, and Exercising Metabolic Rates

Background | Set-up | Lab | Appendix

*Note: The lab presented here is intended for evaluation purposes only. iWorx users should refer to the User Area on [www.iworx.com](http://www.iworx.com) for the most current versions of labs and LabScribe Software.*



**iWorx Systems, Inc.**

[www.iworx.com](http://www.iworx.com)

iWorx Systems, Inc.  
62 Littleworth Road, Dover, New Hampshire 03820  
(T) 800-234-1757 / 603-742-2492 (F) 603-742-2455

LabScribe is a trademark of iWorx Systems, Inc.

©2016 iWorx Systems, Inc.

## Experiment HE-9: Resting, Active, and Exercising Metabolic Rates

### Background

During the day, we expend energy at levels that are related to the activities that we are performing. The rate at which energy is utilized, or the metabolic rate, can be classified into one of three levels: resting, active, and exercising. The energy expended is usually measured indirectly by measuring the rate at which oxygen is consumed during an activity. This method is known as indirect calorimetry.

Humans usually burn 70-80% of the calories that are burned each day at the resting level. These are the calories needed to maintain life. The measurement of the resting metabolic rate (RMR) is performed under the following conditions:

- The subject should not ingest any food during the 12 hours prior to the test.
- The subject should be physically and mentally relaxed.
- The core temperature of the subject should be normal.
- The temperature of the room in which the test is conducted should be comfortable. A temperature that is a few degrees above or below 24°C (75°F) is suitable.

The calories burned during normal daily activities, like walking, working at a computer, and digesting food, are included in the active level of metabolism. To obtain an accurate measurement of the calories expended during an activity, the resting metabolic rate is subtracted from the total metabolic rate recorded during the activity.

During exercise, calories are expended in proportion to the intensity of the exercise. Again, to accurately determine the amount of energy expended during exercise, the resting metabolic rate is subtracted from the total metabolic rate recorded during the exercise.

### Energy Sources

The energy requirements of the body are met with a mixture of energy derived from carbohydrates, fats, and protein. The activity being performed and the stores of carbohydrates and fats available as energy sources determine the proportion of the three macromolecules that are utilized by the body. At rest, a body derives about 40% of its energy from carbohydrates and 60% from fats. As the intensity of activity and the demand for energy increase, a greater proportion of the energy is usually provided by the oxidation of carbohydrates. At the most intense exercise level, all the energy that is required is usually being supplied by carbohydrates.

Protein is usually not a significant source of energy in the body. Unlike fat and carbohydrates in the form of glycogen, the body has no storage deposits of protein. Proteins are important components of tissues, peptide hormones, and enzymes that are continually being broken down and replaced. However, during periods of exercise that are greater than 90 minutes, it is estimated that protein catabolism provides as much as 15% of the energy required. If the body's supply of stored carbohydrates is low from a prior exercise period, protein catabolism could provide as high as 45% of the energy required. Because of its importance in tissues, the utilization of protein as an energy source could cause damage to these tissues. Severe damage to tissues does occur during long-term starvation, when protein is the principal source of energy.

## Indirect Calorimetry

The expenditure of energy can be measured directly using a body calorimeter to measure the amount of heat produced by the body during the burning of calories. Since the number of calories burned and the heat produced during an activity are directly proportional to the amount of oxygen consumed during that activity, the metabolic rate for that activity can be measured indirectly by determining the amount of oxygen consumed during the activity. If the amount of carbon dioxide produced during the activity is also measured, the types and proportions of energy source being utilized during the activity can also be determined.

When examining the oxidation of the three macromolecules used as energy sources, it is shown that the stoichiometries of the reactions for each type of molecules are significantly different.

The oxidation of a molecule of carbohydrate is expressed by the following equation:



As shown in this equation, 6 molecules of carbon dioxide are produced for every 6 molecules of oxygen consumed during the oxidation of carbohydrates, a ratio of 1.0.

The oxidation of a molecule of fatty acid is expressed by this equation:



As shown in this equation, 16 molecules of  $\text{CO}_2$  are produced for every 23 molecules of  $\text{O}_2$  consumed during the oxidation of fatty acids, a ratio of 0.7.

As the ratio of energy supplied by fats and carbohydrates shifts during changes in activity, the ratio of carbon dioxide produced to oxygen consumed also shifts. The ratio of carbon dioxide produced to oxygen consumed during cellular metabolism can be measured by determining the changes in the concentrations of oxygen and carbon dioxide in the air that passes into and out of the lungs. These measurements are possible because the amounts of oxygen and carbon dioxide exchanged between the alveoli and the capillaries in the lungs are directly dependent on the amounts of carbon dioxide produced and oxygen consumed during cellular respiration.

The amounts of oxygen consumed and carbon dioxide produced are measured using an oxygen/carbon dioxide gas analyzer connected to a spirometer. The gas analyzer measures the concentration of oxygen and carbon dioxide in inspired and expired air, and the spirometer determines the volumes of inspired and expired air. When the concentrations and volumes are brought together in a series of equations, the volume of oxygen consumed per minute, known as  $\text{VO}_2$ , and the volume of carbon dioxide produced per minute, known as  $\text{VCO}_2$ , are determined. The ratio of  $\text{VCO}_2/\text{VO}_2$  is the Respiratory Exchange Ratio (RER), which can be used to determine the proportion of carbohydrates and fats utilized, and the energy expended per liter of oxygen consumed, during an activity ([Table HE-9-B1](#)).

The fat and carbohydrate percentages utilized during an activity are determined using the following equations:

$$((1.00 - \text{RER}) / (1.00 - 0.70)) \times 100 = \% \text{ Fat utilized}$$

$$100\% - \% \text{ Fat utilized} = \% \text{ CHO utilized}$$

When heart rates are recorded in conjunction with the rates of oxygen consumption and carbon dioxide production, the effectiveness of different activities to burn calories from fat or carbohydrates, or total calories can be determined by examining the RER at different heart rates. For example, if the RER of a subject is 0.75 at 85 beats per minute (BPM) and 0.90 at 115 BPM, the subject more effectively burns fat at 85 BPM than 115BPM ([Table HE-9-B1](#)).

### ***Energy Expenditure***

The energy expended during an activity is calculated from the RER and the volume of oxygen consumed. For example, if the RER is 0.90, the energy expended is 4.92 kcal/liter O<sub>2</sub>. If 2.5 liters of oxygen are consumed per minute for 20 minutes, a total of 246 kcal are expended during the activity:

$$(2.5 \text{ LO}_2/\text{minute})(20 \text{ min})(4.92 \text{ kcal/liter O}_2) = 246 \text{ kcal}$$

At less intense activity levels, the rates of energy expenditure and RER values are lower. To expend the same amount of energy at a less intense level of activity, the duration of activity must be longer. For example, if the RER is 0.80, the energy expended is 4.80 kcal/liter O<sub>2</sub>. If 1.7 liters of oxygen are consumed per minute, 8.16 kcal are expended per minute:

$$(1.7 \text{ LO}_2/\text{minute})(4.80 \text{ kcal/liter O}_2) = 8.16 \text{ kcal/min}$$

To expend 246 kcal at a rate 8.16 kcal/min would require 30 minutes, 9 seconds:

$$246 \text{ kcal}/(8.16 \text{ kcal/min}) = 30.15 \text{ minutes}$$

**Table HE-9-B1: Respiratory Exchange Ratio (RER) as a Function of the Proportions of Energy Sources.**

RER	Energy (kcal/liter O <sub>2</sub> )	% Energy from CHO	% Energy from Fats
0.70	4.69	0	100
0.75	4.74	15.6	84.4
0.80	4.80	33.4	66.6
0.85	4.86	50.7	49.3
0.90	4.92	67.5	32.5
0.95	4.99	84.0	16.0
1.00	5.05	100	0

In this experiment, students will measure the  $\text{VO}_2$ ,  $\text{VCO}_2$ , and heart rate of a subject while he or she is resting while reclining; performing an activity, like reading or working on a computer while sitting; and, exercising at various levels of intensity from easy to moderate. From these measurements, the subject's RER, proportion of energy sources utilized, and effectiveness of cardiac fitness will be determined.

These measurements will be performed quickly and easily using an iWire-GA  $\text{CO}_2/\text{O}_2$  gas analyzer connected to a spirometry system and a heart rate monitor.

iWorx Sample Lab

## Experiment HE-9: Resting, Active, and Exercising Metabolic Rates

### Equipment Required

PC or Mac Computer

IXTA Data acquisition unit, power supply and USB cable

PHRM-220 Polar<sup>TM</sup> heart rate monitor

FH-1000 Spirometer flow head and plastic tubes

iWire-GA CO<sub>2</sub>/O<sub>2</sub> Gas Analyzer with filter, sample tubing

A-GAK-201 Reusable mask, head gear, and non-rebreathing valve

A-CAL-150 Calibration kit

Flexible tubing, 61 cm long, 3.5cm I.D. (2)

Rubber adapter, 3.5cm O.D. to 2.85cm I.D.

5 Liter mixing chamber

Treadmill or exercise bike

3 Liter Calibration syringe

### Setup the IXTA and iWire-GA

1. Connect the iWire-GA to the iWire1 port on the front of the IXTA, and plug it into the wall using the power supply.
2. Plug the IXTA into the wall and, using the USB cable, to the computer.

**NOTE: The iWire-GA must be plugged into the IXTA prior to turning both machines on.**

3. Turn on the IXTA and the iWire-GA.
4. Open LabScribe.
5. Click Settings → Human Exercise-iWireGA → RestActExerMetabolism.
6. Once the settings file has been loaded, click the **Experiment** button on the toolbar to open any of the following documents:
  - Appendix
  - Background
  - Labs
  - Setup (opens automatically)

### Setup the Metabolic Cart

1. Locate the A-FH-1000 flow head and tubing in the iWorx kit ([Figure HE-9-S1](#)).



Figure HE-9-S1: The A-FH-1000 flow head, and airflow tubing.

2. Carefully attach the two airflow tubes onto the two sampling outlets of the A-FH-1000 flow head and the other ends of the two airflow tubes onto Channel A1 on the front of the IXTA ([Figure HE-9-S4](#)).

**Note:** Make sure to connect the airflow tubing so that the ribbed tube is attached to the red outlet port of the flow head and also to the red inlet port of the spirometer. The smooth side of the tubing attaches to the white ports.

3. Locate the mixing chamber in the iWorx kit ([Figure HE-9-S2](#)).
4. Connect the inlet of the A-FH-1000 flow head to the outlet of the mixing chamber ([Figure HE-9-S3](#)).

**Note:** Be sure to connect the flow head to the mixing chamber so that the red outlet port is facing towards the mixing chamber.





Figure HE-9-S2 and HE-9-S3: The mixing chamber showing the 1000L/min flow head connected to the outlet.




5. Locate the non-rebreathing valve, mask, and smooth interior tubing in the iWorx kit ([Figure HE-9-S5](#)).
6. Attach one end of the smooth interior tubing to the inlet of the mixing chamber ([Figure HE-9-S6](#)), and the other end to the outlet of the non-rebreathing valve. There are arrows on the valve that indicate the direction of air flow.
7. Attach the mask to the side port of the non-rebreathing valve.



Figure HE-9-S4: The iWire-GA gas analyzer connected to an IXTA. All tubings are connect properly in this image.





*Figure HE-9-S5: Mask, non-rebreathing valve, and smooth interior tubing.*

8. On the iWire-GA, place one filter on the “Room Air” port, place a second filter on the “Sample In” port. Attach the braided end of the Nafion sampling tube to the filter on the “Sample In” port.
9. Place the other end of the Nafion sampling tube on the gas sampling port near the outlet of the mixing chamber ([Figure HE-9-S6](#)).



*Figure HE-9-S6: The assembled devices used during metabolic studies. The assembly includes: the mixing chamber, smooth interior tubing, Nafion sampling tubing, flow head, spirometer, non-rebreathing valve, and mask.*

10. Plug the outlet tubing from the iWire-GA to the port on the mixing chamber, opposite the flow head.
11. Locate the PHRM-220 Polar<sup>™</sup> heart rate monitor transmitter, electrode belt, and receiver in the iWorx kit ([Figure HE-9-S8](#)).

12. Plug the connector of the PHRM-220 receiver into the “HR” input on the iWire-GA.



*Figure HE-9-S8: The PHRM-220 transmitter, belt, and receiver set.*

**The non-rebreathing valve can be used with the attached mask or with an optional mouthpiece.**

**If the subject is using a mask (preferred method):**

- Attach the head gear to the mask.
- Attach the non-rebreathing valve to the mask. Depending on the model of the mask, an adapter may be required.
- Instruct the subject to try on the assembly. Adjust the straps so that the mask fits the subject comfortably. Make sure there are no leaks around the mask.

**If the subject is using a mouthpiece:**

- Attach the headgear to the brackets on the non-rebreathing valve. The pair of straps with the narrowest spacing go over the top of the subject's head.
- Connect the mouthpiece to the side port of the valve so that the valve is oriented horizontally, and the saliva trap of the mouthpiece is pointed downward.
- Instruct the subject to try on the assembly. Adjust the straps so that the mouthpiece fits the subject comfortably. Make sure there are no leaks between the mouthpiece and the valve or around the mouthpiece.

**The gas analyzer must warm up for at least 15 minutes.**

***Note: For increased accuracy, users must complete the flow Head Calibration procedure. Please see Appendix I for directions on how to perform this calibration. The calibration of the 1000L flow head requires a 3L Calibration Syringe.***

## Load a PreSaved flow Head Calibration (\*.iwxgcd) File

**Note:** This procedure is used once a calibration curve has been generated using the Spirometer Calibration directions in Appendix I.

All of the following directions will be prompted by the software. Follow the directions as they pop up on the LabScribe software.

1. Load the lab settings file you wish to perform as stated in the “Start the Software” section ([Figure HE-9-S9](#)).
2. Assemble the spirometer, flow head, tubing, mixing chamber and calibration syringe as shown in Appendix I or in the SpirometerCalibration directions.
3. Click the Setup button shown in the left side window. Follow the directions as prompted by the ([Figure HE-9-S10](#)). The Online Setup Dialog window will open
  - Enter your subject's information or Load a subject from a previously saved file.
  - Click “Settings” to change any parameters you wish to view ([Figure HE-9-S11](#)).




Figure HE-9-S9: Initial screen for starting a Fitness Assessment test. Follow the directions as prompted by the buttons on the left side of the window.

Online Metabolic Setup Dialog

Subject Settings

Load Subject

Name

Protocol

Age

Sex

Height(cm)

Weight(kg)

Maximum Heart Rate

Blood Pressure  /

OK Cancel

Figure HE-9-S10: Online Metabolic Setup Dialog window.

Online Metabolic Setup Dialog

Subject Settings

CO2 Channel

O2 Channel

Volume Channel

Average  seconds

O2 Conc. in Inhaled Air  %

CO2 Conc. in Inhaled Air  %

Calculations

- Sel
- Time
- Abs.VO2
- Abs.VCO2
- Rel.VO2
- Rel.VCO2
- RER
- REE
- METS
- %Fat
- %Carbs
- Calories
- Min O2 conc
- Max CO2 conc
- Expired Vol.

OK Cancel

Figure HE-9-S11: Settings dialog of the Online Metabolic Setup window.

4. Perform the Quick Flow Calibration by clicking the button and following the prompted directions ([Figure HE-9-S12](#)).
  - When you click “Load”, you will be prompted to load the .iwxcd file created when you performed the full flow head calibration.




Figure HE-9-S12: Spirometer Calibration dialog window.

5. Select Save As in the File menu, type a name for the file. Click on the Save to save the data file.

### Calibrating the iWire-GA Gas Analyzer

**Note:** Warm up the gas analyzer for at least 15 minutes prior to use. Make sure the calibration gas tank is located close to the gas analyzer.

This procedure will calibrate both the O<sub>2</sub> and CO<sub>2</sub> channels.

Connect the gas sample tubing of the A-CAL-150 Calibration Kit ([Figure HE-9-S13](#)) to the Luer-Lock connector on the output of the regulator.



Figure HE-9-S13 Calibration Kit (A-CAL-150).



1. Click the Calibrate Gas Analyzer button. Click Perform Quick Software Gas Calibration.
2. Follow the directions as prompted. Room air will be sampled for 10 seconds. Calibration gas will be sampled for 15 seconds.
3. If necessary, move the cursors into correct position ([Figure HE-9-S14](#)).




Figure HE-9-S14: Advanced units conversion dialog for room air and calibration gas.



## Appendix I: Initial Spirometer flow head Calibration


**For accuracy of measurements, users must include this calibration procedure as part of the Exercise Physiology Lab protocol.**

***It is suggested that this procedure be followed at the beginning of every term and when using a new flow head-spirometer combination.***


***Note:*** This calibration protocol precedes the actual calibration of the GA-200 or GA-300 gas analyzer. You will not need the gas analyzer at this time.

***Note:*** Whenever you will be using a different flow head, you will need to repeat this calibration procedure from the beginning by loading a new Spirometer Calibration settings file.

1. Open the LabScribe software.
2. Click Settings - Human Exercise-iWireGA. Choose SpirometerCalibration to launch the calibration settings file.
3. Assemble the flow head, tubing, mixing chamber and calibration syringe.
4. Plug the tubing into the internal spirometer channel A1.
  - Connect the flow head to the IXTA using the flow head tubing, making sure that the ribbed side of the tubing connects the red marked port on the flow head to the red marked port on the spirometer ([Figure HE-9-S15](#)).
  - Connect the smooth side of the tubing to the other ports.
5. Connect end of the 1000L flow head with the red marked onto white flange of the mixing chamber. Make sure the tubing is in an upright direction ([Figure HE-9-S16](#)).




*Figure HE-9-S15: The 1000L flow head and Figure HE-9-S16: The 1000L flow head attached to the mixing chamber showing the tubing in an upright position and the red port facing the mixing chamber.*



*Figure HE-9-S17: 1000L flow head connected to the mixing chamber, showing the Nafion tubing connected to the outlet sampling port near the flow head.*

**Note:** Make sure the red port on the flow head faces into the mixing chamber.

6. Connect one end of the smooth bore tubing to the 3L calibration syringe as shown in [Figure HE-9-S18](#).
7. Connect the other end of the smooth bore tubing to the mixing chamber, opposite the flow head ([Figure HE-9-S19](#)).



*Figure HE-9-S18: 3 liter calibration syringe connected to the smooth bore tubing and Figure HE-9-S19: The smooth bore tubing connected to the mixing chamber.*

8. If also setting up the gas analyzer at this time:
  - Connect the braided Nafion tubing to the filter on the gas analyzer and to the flow head side of the mixing chamber. Make sure the braided end is connected to the filter ([Figure HE-9-S17](#)).
  - Connect the thin flexible tubing from the outlet of the gas analyzer to the port next to the smooth bore tubing on the opposite side of the mixing chamber.

9. If not using the gas analyzer at this time, connect the flexible tubing from the port on one side of the mixing chamber to the port on the other. This ensures there is no air leaking from the chamber.
10. Pull the plunger on the 3L Calibration Syringe all the way out until it stops.
11. Click the Record button.
12. Wait for at least 10 seconds of recording so that there is no flow of air moving through the syringe.
13. Push the plunger in all the way until it stops. Pull the plunger out all the way until it stops.
14. Repeat the procedure in Step 13, for at least 50 repetitions, varying the speed and force on the plunger. Make sure to pause between strokes.
15. The faster the speed of the stroke, the higher the flow through the calibration syringe.

**Note:** Ideally the flow head calibration recording should span air flow values to include the minimum to maximum flow levels for the particular experiment being conducted.

16. After at least 50 repetitions have been performed, wait at least 5 seconds after the final repetition and then click Stop.
17. Select Save As in the File menu, type a name for the file.
18. Click on the Save button to save the raw data for generation of a flow head calibration \*.iwxgcd file.
19. Click AutoScale on the Air flow channel.
20. Use the Display Time icons to adjust the Display Time of the Main window to show the complete calibration data ([Figure HE-9-S20](#)).
21. Click the Double Cursor icon so that two cursors appear on the Main window.




Figure HE-9-S20: The LabScribe toolbar.

22. Click Advanced on the main toolbar. Then click Metabolic, and Calibrate flow head ([Figure HE-9-S21](#)).
23. Place the two blue vertical cursors so that:

- The left-hand most cursor is on the flat line prior to the start of the calibration data. Make sure the cursor is at the beginning of the 10 second baseline.
- The right-hand most cursor is on the flat line after the final calibration stroke ([Figure HE-9-S22](#)).




Figure HE-9-S22: Calibrate flow head dialog window.

24. In the new window that opens ([Figure HE-9-S24](#)), enter these values:

- flow channel = Expired Air flow
- Baseline = Use the first 10 seconds as zero
- Calibrate difference between cursors to 3 L.

25. Click the Calibrate the difference between cursors to button. This will generate the curve as shown above.

26. A new window will open prompting you to Save your file as an \*.iwxgcd flow head calibration file. Name your file and click Save.

27. Click OK.

**Note:** At this point, a raw calibration data file (\*.iwxdata) and a flow head calibration file (\*.iwxgcd) have been generated.

28. Exit LabScribe or open a Human Exercise lab settings file.

**Note:** Once a saved \*.iwxgcd file is loaded, a simple 5-10 stroke calibration procedure can be used to update the file for immediate use.






Figure HE-9-S23 The calibration recording showing the vertical cursors in the correct position for generating a calibration curve. Note – the recording you generate should look similar to this.




Figure HE-9-S24: Calibration syringe data.

## Experiment HE-9: Resting, Active, and Exercising Metabolic Rates

### Before Starting

1. Read the procedures for the experiment completely before beginning the experiment. Have a good understanding of how to perform the experiment before making recordings.
2. It is important that the subject is healthy and has no history of respiratory or cardiovascular problems.
3. Allow the spirometer to warm up for 15 minutes before recording for the first time.
4. Determine if the airflow tubes between the flow head to the spirometer amplifier are attached to the proper inlets on each device.
  - Since this test does not need to be recorded, click on the Save to Disk button in the lower left corner of the Main window. If LabScribe is in Preview mode, there will be a red X across the Save to Disk button.
  - Click on the Preview button.

**Note:** If the user clicks the Preview button and an error window appears on the Main window indicating the iWorx hardware cannot be found, make sure the iWorx unit is turned on and connected to the USB port of the computer. Then, click on the OK button in the error window. Pull down the LabScribe Tools menu, select the Find Hardware function, and follow the directions on the Find Hardware dialogue window.

- Have the subject inhale and exhale through the mask 2 or 3 times while the complete spirometry circuit is assembled.
  - Click on the AutoScale button at the upper margin of the Expired Air Flow and Lung Volume channels.
  - If the proper end of the flow head is attached to the outlet of the mixing chamber, the traces on the Air Flow and Lung Volume channels will go up when the subject exhales.
  - If the traces on these channels go down during exhalation, remove the flow head from the outlet of the mixing chamber and place the other end of the flow head on the outlet of the mixing chamber.
  - Click on the Stop button.
5. Click on the Save to Disk button, in the lower left corner of the Main window, to change LabScribe from Preview mode to Record mode. If LabScribe is in Record mode, there will be a green arrow on the Save to Disk button.



## Exercise 1: Resting Metabolic and Heart Rates

Aim: To determine the amount of oxygen consumed ( $\text{VO}_2$ ), carbon dioxide produced ( $\text{VCO}_2$ ), respiratory exchange ratio (RER) in relation to the heart rate of a resting subject.

### Procedure

1. Instruct the subject to sit quietly, become accustomed to breathing through the spirometry equipment, and breathe normally before any recordings are made.
2. Once the subject and recording equipment are all prepared, disconnect the smooth-bore tubing from the mixing chamber to ensure that no air is entering the system at this time.
3. Type <Subject's Name> Inhalation at Rest in the Mark box that is to the right of the Mark button.
4. Click on the Record button. After waiting ten seconds for the Expired Air Volume ATP channel to zero, reconnect the proper end of the tubing to the mixing chamber. Press the Enter key on the keyboard to mark the recording.

*Note: The LabScribe software will zero the Expired Air Volume ATP channel during the first ten seconds of recording. No air should be moving through the flow head during this time.*

5. Click the AutoScale buttons on all channels.
6. On the Expired  $\text{CO}_2$  Concentration (%) channel, notice that the  $\text{CO}_2$  concentration increases in the first few minutes of the recording and then reaches a near-steady level.
  - The time that it takes the chamber to be filled with expired air and reach a near-steady level of carbon dioxide is dependent on the tidal volume and respiration rate of the subject and the volume of the mixing chamber. It will take longer to fill the chamber if the subject's respiration rate and tidal volume are low, or the chamber is large.
  - Every breath exhaled into the mixing chamber pushes a matching volume of expired air out of the mixing chamber.
  - Record baseline data, while the mixing chamber air is replaced with the subject's expired air, for approximately 5-10 minutes prior to beginning any experiments.
7. On the Expired  $\text{O}_2$  Concentration (%) channel, notice that the  $\text{O}_2$  concentration decreases in the first few minutes of the recording and then stays a near-steady level. As pointed out in the previous step, the size of the mixing chamber, the tidal volume, and respiration rate of the subject, determine the time it takes for the concentration of oxygen in expired air to reach that near-steady level.
8. On the Expired Air Volume channel, the STPD Volume-MC spirometry function converts the data from the Expired Air Flow channel to the volumes of expired air at STPD. Notice that the recorded volume increases in a ramp-like manner with each breath.




Figure HE-9-L2: Gas concentrations and volumes of a standing subject displayed on the Main window. Concentrations of gases reach a steady state after the mixing chamber is filled. Cursors placed one minute apart indicate  $V\text{CO}_2$  is 0.284 liters/min and  $\text{VO}_2$  is 0.397 liters/min. Therefore, RER is 0.96.

9. Continue to record until one or more minutes of data are recorded while the concentrations of oxygen and carbon dioxide in expired air are at a steady level. Once the appropriate duration of data is recorded, click Stop to halt the recording. Your recording should be similar to the data displayed in [Figure HE-9-L2](#).
10. Select Save As in the File menu, type a name for the file. Choose a destination on the computer in which to save the file, like your lab group folder). Designate the file type as \*.iwxdata. Click on the Save button to save the data file.

### Data Analysis

1. Display the complete data recording in the Main window. Use the Display Time icons to adjust the Display Time of the Main window to show the complete recording on the Main window.

**Note:** To simplify the display of data on the Main and Analysis windows, select an alternate view of the window. Alternate views allows channels to be hidden from view on a window even though they may still be used in computations that create channels still in view.

- Open the View menu, which is located to the right of the 1-Cursor icon on the LabScribe toolbar, by clicking on the arrow to the right of the current view, Default View.
  - Select VO<sub>2</sub>-VCO<sub>2</sub> from the list on the menu. This view is programmed in the settings file for this experiment. The Expired Air Volume ATP and Expired Air Volume STPD channels are hidden from view, but still used in the calculations for VO<sub>2</sub> and VCO<sub>2</sub>.
2. Select and display the last 60-second section of the recording while the oxygen and carbon dioxide concentrations were at a steady level on the Main window. Select the 60-second section of the recording by:
    - Placing the cursors on either side of the 60-second section of data; and
    - Clicking the Zoom between Cursors button on the LabScribe toolbar ([Figure HE-9-L3](#)) to expand the selected section of data to the width of the Main window.
  3. Click on the Analysis window icon in the toolbar or select Analysis from the Windows menu to transfer the data displayed in the Main window to the Analysis window ([Figure HE-9-L4](#)).




Figure HE-9-L3: The LabScribe toolbar.

4. Look at the Function Table that is above the uppermost channel displayed in the Analysis window. The mathematical functions, V<sub>2</sub>-V<sub>1</sub>, T<sub>2</sub>-T<sub>1</sub>, and Mean, should appear in this table. Values for these two parameters on each channel are seen in the table across the top margin of each channel.
5. Once the cursors are placed in the correct positions for determining the carbon dioxide and oxygen concentrations in expired air, the expired air minute volume at STPD, the VCO<sub>2</sub>, and the VO<sub>2</sub> of the resting subject, the values for these parameters can be recorded in the on-line notebook of LabScribe by typing their names and values directly into the Journal.
6. The functions in the channel pull-down menus of the Analysis window can also be used to enter the names and values of these parameters from the recording to the Journal. To use these functions:
  - Place the cursors at the locations used to measure the concentrations and volumes.
  - Transfer the names of the mathematical functions used to determine these concentrations and volumes to the Journal using the Add Title to Journal function in the VCO<sub>2</sub> Channel pull-down menu.

- Transfer the values for the volumes and rates to the Journal using the Add All Data to Journal function in the VCO<sub>2</sub> Channel pull-down menu.
7. Use the mouse to click on and drag the cursors to positions on the Expired O<sub>2</sub> Concentration (%) channel that are one minute apart during the near-steady level of oxygen in the mixing chamber ([Figure HE-9-L4](#)). The values for the following parameters are determined when the cursors are positioned as directed:
    - Mean concentration of CO<sub>2</sub> in expired air, which is the value for Mean on the Expired CO<sub>2</sub> Concentration channel.
    - Mean concentration of O<sub>2</sub> in expired air, which is the value for Mean on the Expired O<sub>2</sub> Concentration channel.
    - Air Volume at STPD expired in a minute, which is the value for V2-V1 on the Expired Air Volume STPD channel.
    - Volume of carbon dioxide produced in one minute (VCO<sub>2</sub>), which is the value for V2-V1 on the VCO<sub>2</sub> channel.
    - Volume of oxygen consumed in one minute (VO<sub>2</sub>), which is the value for V2-V1 on the VO<sub>2</sub> channel.
    - Mean heart rate during the period, which is the value for Mean on the Heart Rate channel.
  9. Record the values in the Journal using one of the techniques described in Steps 6 or 7.
  10. Record the values for the mean CO<sub>2</sub> and O<sub>2</sub> concentrations in expired air, the minute volume of expired air at STPD, the VCO<sub>2</sub>, the VO<sub>2</sub>, and the heart rate while the subject was resting in [Table HE-9-L1](#).
  11. Calculate the Respiratory Exchange Ratio (RER) of the resting subject by dividing the value for VCO<sub>2</sub> by the value for VO<sub>2</sub>. Record the values for RER in the data table.
  12. Calculate the Relative VO<sub>2</sub> of the resting subject by dividing the value for VO<sub>2</sub> by the subject's body weight (kg). Record the values for Relative VO<sub>2</sub> in the data table.
  13. Calculate the VO<sub>2</sub>/Heart Rate Ratio (VO<sub>2</sub>/HR) of the resting subject by dividing the value for Absolute VO<sub>2</sub> by the heart rate (HR) for the same period. Record the value for VO<sub>2</sub>/HR in the data table.




Figure HE-9-1A: Gas concentrations, air volumes, VCO<sub>2</sub>, VO<sub>2</sub>, and heart rate of a resting subject displayed on the Analysis window. Cursors placed one minute apart indicate VCO<sub>2</sub> is 0.284 liters/min, the VO<sub>2</sub> is 0.297 liters/min, and the mean heart rate is 72.044. Therefore, RER is 0.96 and the VO<sub>2</sub>/HR is 4.122ml/min/beat.

## Exercise 2: Active Metabolic and Heart Rates

Aim: To determine the amount of oxygen consumed (VO<sub>2</sub>), carbon dioxide produced (VCO<sub>2</sub>), respiratory exchange ratio (RER) in relation to the heart rate of a active subject.

### Procedure

1. Use the same procedures used in Exercise 1 to record the lung volumes from the subject while he or she is performing a typical daily activity.
2. Remember to record until at least two or more minutes of data are recorded while the concentrations of oxygen and carbon dioxide in expired air are at a steady level. Once the appropriate duration of data is recorded, click Stop to halt the recording.
3. Mark the recording with comments that indicate the name of the subject and the beginning of each minute being recorded.



### **Data Analysis**

1. Use the same procedures used in Exercise 1 to determine the minute volumes of expired air, oxygen consumed ( $\text{VO}_2$ ), and carbon dioxide produced ( $\text{VCO}_2$ ) during the activity.
2. Record the values for the mean  $\text{CO}_2$  and  $\text{O}_2$  concentrations in expired air, the minute volume of expired air at STPD, the  $\text{VCO}_2$ , the  $\text{VO}_2$ , and the heart rate while the subject was resting in [Table HE-9-L1](#).
3. Calculate the Respiratory Exchange Ratio (RER) of the resting subject by dividing the value for  $\text{VCO}_2$  by the value for  $\text{VO}_2$ . Record the values for RER in the data table.
4. Calculate the Relative  $\text{VO}_2$  of the resting subject by dividing the value for  $\text{VO}_2$  by the subject's body weight (kg). Record the values for Relative  $\text{VO}_2$  in the data table.
5. Calculate the  $\text{VO}_2$ /Heart Rate Ratio ( $\text{VO}_2/\text{HR}$ ) of the resting subject by dividing the value for Absolute  $\text{VO}_2$  by the heart rate (HR) for the same period. Record the value for  $\text{VO}_2/\text{HR}$  in the data table.

### **Exercise 3: Exercising Metabolic and Heart Rates**

**Aim:** To determine the amount of oxygen consumed ( $\text{VO}_2$ ), carbon dioxide produced ( $\text{VCO}_2$ ), respiratory exchange ratio (RER) in relation to the heart rate of a subject performing light to moderate exercise.

#### **Procedure**

1. Use the same procedures used in Exercise 1 to record the lung volumes from the subject while he or she is exercising at an easy level, like 3 km/hour (2 miles/hr) on a treadmill or 10 km/hour (6 miles/hr) on an exercise bike.
2. Remember to record until at least two or more minutes of data are recorded while the concentrations of oxygen and carbon dioxide in expired air are at a steady level.
3. Continue to record as the exercise level is doubled to a higher level, like 6 km/hour (4 miles/hr) on a treadmill or 20 km/hour (12 miles/hr) on an exercise bike. Remember to record until one or more minutes of data are recorded at a steady level.
4. Continue to record as the exercise level is increased to a higher level, like 9 km/hour (6 miles/hr) on a treadmill or 30 km/hour (18 miles/hr) on an exercise bike. Remember to record until one or more minutes of data are recorded at a steady level.

**Warning:** *The subject should stop exercising if they begin to feel any discomfort.*

5. Mark the recording with comments that indicate the name of the subject, the level of exercise being performed, and the beginning of exercise level being recorded.



**Table HE-9-L1:  $\text{VCO}_2$ ,  $\text{VO}_2$  and the Respiratory Exchange Ratio (RER) of a Subject at Resting, Active, and Exercising Levels of Metabolism.**

Room Temperature_____ °C	Activities					
	Units	Resting	Daily Activity	Exercise Level 1	Exercise Level 2	Exercise Level 3
Barometric Pressure_____ mmHg						
Exhaled [Carbon Dioxide]	% $\text{CO}_2$					
Exhaled [Oxygen]	% $\text{O}_2$					
Expired Minute Volume STPD	L/min					
Absolute $\text{VCO}_2$	ml/min					
Absolute $\text{VO}_2$	ml/min					
Heart Rate	beats/min					
RER						
Relative $\text{VO}_2$	ml/min/kg					
$\text{VO}_2/\text{HR}$	ml/beat					

### Data Analysis

1. Use the same procedures used in Exercise 1 to determine the minute volumes of expired air, oxygen consumed ( $\text{VO}_2$ ), and carbon dioxide produced ( $\text{VCO}_2$ ) during the activity.
2. Record the values for the mean  $\text{CO}_2$  and  $\text{O}_2$  concentrations in expired air, the minute volume of expired air at STPD, the  $\text{VCO}_2$ , the  $\text{VO}_2$ , and the heart rate while the subject was resting in the above data table.
3. Calculate the Respiratory Exchange Ratio (RER) of the resting subject by dividing the value for  $\text{VCO}_2$  by the value for  $\text{VO}_2$ . Record the values for RER in the above data table.
4. Calculate the Relative  $\text{VO}_2$  of the resting subject by dividing the value for  $\text{VO}_2$  by the subject's body weight (kg). Record the values for Relative  $\text{VO}_2$  in the data table.
5. Calculate the  $\text{VO}_2/\text{Heart Rate Ratio}$  ( $\text{VO}_2/\text{HR}$ ) of the resting subject by dividing the value for Absolute  $\text{VO}_2$  by the heart rate (HR) for the same period. Record the value for  $\text{VO}_2/\text{HR}$ .

## Questions

1. During which activity was the subject's RER the highest? In which period was it the lowest?
2. During which activity was the subject consuming a higher proportion of fat calories?
3. During which activity was the subject consuming a higher proportion of carbohydrate calories?
4. During which activity did the subject consume the highest amount of total calories for the level of cardiac activity?
5. During which activity did the subject consume the highest amount of fat calories for the level of cardiac activity?
6. Evaluate the physical fitness of your subject. How does the level of your subject's physical fitness correlate to his or her RER while resting, being active, and exercising?
7. Evaluate the diet of your subject. How does your subject's diet correlate to his or her RER while resting, being active, and exercising?
8. How does your subject's level of physical fitness, diet, RER, and cardiac effectiveness correlate to those parameters from other members of the class?

## Calculations: Calculate and Plot Metabolic Parameters

Values for  $\text{VO}_2$ ,  $\text{VCO}_2$ , RER, TV, and other parameters ([Table HE-9-L2](#)) from the segments of the test can be calculated automatically by using the Metabolic Calculations window.

1. To use the Metabolic Calculations window, pull down the Advanced menu and select Metabolic. Select Mixing Chamber: Offline Calculations from the submenu to open the Metabolic Calculations Dialog window.
2. On the left side of the Metabolic Calculations window:
  - Pull down the  $\text{CO}_2$ ,  $\text{O}_2$ , Volume, Heart Rate, and Energy Channel menus to select the channels on which the  $\text{CO}_2$  and  $\text{O}_2$  concentrations, lung volumes, heart rates, and workload were recorded.
  - When analyzed, the data file will be divided into time segments. The average of each parameter in each segment will be reported in the data table on the Metabolic Calculations window. Enter the time (in secs) in the Average box to select the time length of each segment.
  - In the  $\text{O}_2$  and  $\text{CO}_2$  Concentrations in Inhaled Air boxes, enter the concentrations of oxygen and carbon dioxide in the inhaled air, which is room air in most tests.
4. Click on the Calculate button on the left side of the Metabolic Calculations Dialog window to calculate the average value of each parameter listed in the table for each time segment of the recorded data, and to plot the selected parameters against each other in the plot panel ([Figure HE-9-L5](#)).
5. In the lower left corner of the plot panel, click on the arrow to open the pull-down menu listing the types of plots ([Table HE-9-L3](#)) that can be made with the metabolic parameters calculated by this analytical tool. Select the plot to be displayed in the plot panel when the calculations are performed.

**Table HE-9-L2: List of Parameters Calculated on the Mixing Chamber Offline Metabolic Window**

Term	Parameter	Description	Units
Abs.VO <sub>2</sub>	Absolute VO <sub>2</sub>	Volume of oxygen (O <sub>2</sub> ) consumed / minute	Liters/minute
Abs.VCO <sub>2</sub>	Absolute VCO <sub>2</sub>	Volume of carbon dioxide (CO <sub>2</sub> ) produced per minute	Liters/minute
Rel.VO <sub>2</sub>	Relative VO <sub>2</sub>	Volume of O <sub>2</sub> consumed per kg body weight per minute	ml/kg/minute
Rel.VCO <sub>2</sub>	Relative VCO <sub>2</sub>	Volume of CO <sub>2</sub> produced per kg body weight per minute	ml/kg/minute
RER	Respiratory Exchange Ratio	Ratio of VCO <sub>2</sub> /VO <sub>2</sub>	None
REE	Resting Energy Expenditure	5.46 (AbsoluteVO <sub>2</sub> ) + 1.75 (AbsoluteVCO <sub>2</sub> )	kcal/day
METS	Metabolic Equivalent of Task	1 MET = 3.5ml O <sub>2</sub> /kg/min or 1kcal/kg/hr	MET
O <sub>2</sub> Min.	O2 Minimum - exhalation	Minimum concentration of O <sub>2</sub> recorded during test period	Percentage
CO <sub>2</sub> Max.	CO2 Maximum - exhalation	Maximum concentration of CO <sub>2</sub> recorded during test period	Percentage
VE	Expired Tidal Volume	Volume of air displaced during normal exhalation	Liters/breath
P	Power	Workload during the stages of the test	Watts
HR	Heart Rate	Number of beats in a minute: (60 sec/min) divided by the (secs/breath)	Beats / Minute

**Table HE-9-L3: Plots Available on the Offline Metabolic Window.**

	Available Plots										
Y-Axis Parameter 1	VO <sub>2</sub>	VCO <sub>2</sub>	V <sub>e</sub>	V <sub>e</sub>	HR	V <sub>t</sub>	V <sub>e</sub>	HR	VO <sub>2</sub>	V <sub>e</sub> /VO <sub>2</sub>	RER
Y-Axis Parameter 2	VCO <sub>2</sub>				VCO <sub>2</sub>			VO <sub>2</sub> / HR	VCO <sub>2</sub>	V <sub>e</sub> / VCO <sub>2</sub>	
Y-Axis Parameter 3	RER										
X-Axis Parameter	Time	VO <sub>2</sub>	VO <sub>2</sub>	VCO <sub>2</sub>	VO <sub>2</sub>	V <sub>e</sub>	Watts	Watts	Watts	Watts	Watts




Figure HE-9-L5: The metabolic parameters, and plots of  $\dot{V}O_2$ ,  $\dot{V}CO_2$ , and RER vs. Time, displayed in the Metabolic Calculations window used offline to analyze data collected during an aerobic fitness test. Notice that the  $\dot{V}O_2$  and  $\dot{V}CO_2$  values increase quickly as the subject performs more strenuous segments of the test.