

# iWorx ROAM-ECG

## Wireless and Wearable ECG System

The iWorx ROAM-ECG is a versatile, wireless system designed to provide long-term recording and unrestricted movement. The recording device easily attaches to the subject's chest with standard electrodes. The protective case doubles as both a charging station and a data hub. The system works with both Windows and OSX via standard USB interface. Data can be quickly analyzed with the included LabScribe analysis software or exported directly for other processing. The ROAM-ECG allows for hours of continuous recording with an untethered reach of 20 feet.


# ROAM



Over 20Ft



Of Range



## Features

- 8+ Hours of Continuous Recording
- Plug and Play Technology
- PC and Mac
- USB-Powered
- Fully-Modifiable Filters
- Analysis, Data Processing, and Experiment Builder
- No Syncing Required

## Includes

- ROAM-ECG Recorder
- ROAM-Wireless Receiver
- Retractable-USB Cord
- Hard-Shell Case for Docking and Storage
- LabScribe 4.0 Analysis Software

