

Making Physiology Happen

The IX-GSR recording module allows recording of GSR (Galvanic Skin Response) or EDA (Electrodermal Activity). The recorder includes iWorx LabScribe™ Recording and Analysis Software and connects directly to a USB port on either a Macintosh or Windows computer. Attaches to the subject via finger electrodes.

iWorx Isolated GSR Recorder



IX-GSR Specifications:

- Analog Input: 1 Galvanic Skin Response
- GSR Range 0-10,000 micro Siemens

A/D Converter:

- Sampling Speed: 20k Samples/sec per channel simultaneous
- Resolution: 16 Bit
- Interface: USB full speed

General:

- Enclosure: Plastic
- Certification: CE
- Power: USB
- Warranty: Protected by a 3 year warranty
- Software: LabScribe™

Add-Ons:



- TM-BIO Temperature Sensor